



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motor Driven Power Supply

**Product Code :**  
LBNY-0005-18400010



**Description :**

Suitable for Power Supply with variable voltage the braking systems and the excitation of the machines through manual or automatic operation.

Automatic regulation of excitation to keep a constant voltage

DC output: 0 to 210 V, 2 A

Power supply: 220 V, 50/60 Hz

**Technical Specification :**

Motor Driven Power Supply



**Laboratory Instrument India**