




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Maximum Demand Meter

**Product Code :**  
LBNY-0005-18400013



**Description :**

Measurement of voltages, currents, frequencies, active power, reactive power, apparent power. Microprocessor controlled three phase power analyzer.

Input current: 5 A (max 20 Arms)

Input voltage: 500 V (max 800 Vrms)

Auxiliary supply: single phase from mains

Operating frequency: 47 - 63 Hz.

**Technical Specification :**

Maximum Demand Meter



**Laboratory Instrument India**