



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Laboratory Filtration Pump

**Product Code :**  
LBNY-0001-800015



**Description :**

Laboratory Filtration Pump

**Technical Specification :**

Laboratory filtration pump

Laboratory filtration pump is a vacuum pump is a liquid transfer pump. It has a compact design that makes the pump ideal for use on both the workbench and in the laminar flow hood.

General Features:

Compact design which allow for use on both the work bench and in the laminar flow hood.

Allow direct transfer (no handling of liquid waste).

High-performance flow rate

Noise reduction design

Maintenance-free technology

Allow easy decontamination after use and compatible with different decontaminants

Vacuum

Up to 700 mbar

Flow rate (measured with water at 20°C)

3 Liter/min

Temperature (5-40?) and relative humidity (up to 80%)

Altitude:- 2300-2400 meter

Voltage: 220VAC, + 10% , 50Hz

Protection type and class :- IP64, III.



**Laboratory Instrument India**