



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Needle Smelter

**Product Code :**  
LBNY-0001-800016



**Description :**

Needle Smelter

**Technical Specification :**

Needle Smelter

It melts used needles safely.

Needle Smelter as a right measure to dispose used needles to prevent potential needle sticking accidents.

Needle Smelter disposes used needles completely and easily by melting needles with heat of 1,500° to be melted into powder. It can process various syringe needles from 15-31 gauge with any length within one second.

Needle Smelter is small and easy to operate.

Technical Specification:

Capacity: 15 - 31 gauge needles

Disposable time: 80 -100 needles per minute

Durability: Semi-permanent

Weight: Size (W x H x D): (230-250) mm, (100 - 120) mm, (110 - 130) mm

Heating capacity: 1500 oC

Ambient temperature 5 to 40°C

230 V +-10% , /50 Hz

---

# Laboratory Instrument India

