



Email : sales@laboratoryinstrumentindia.com

Product Name :
Antenna Trainer

Product Code :
LBNY-0005-1900003



Description :

The trainer includes: basic panel with generator, Lecherline, RF detector, mobile receiver, supporting base for receiver, reflector, joint for fixing to the base, locking knob, shielded connecting cable between the Lecher line and the detector with banana plugs, diam. 2 mm (length = 60cm), coaxial support for antenna(mast) with turning base, set of cables for measuring and connections, T BNC, different lengths, ground plane antenna, Yagi antenna, set of 5 monopoles with banana terminations, coaxial cable RG59 with F/F BNC, folded dipole antenna, simple dipole antenna, length = 60 cm complete with M/M and F/F adapters.

Input power supply: - 15 Vdc, 200 mA

Frequency: from 860 to 950 MHz.

Technical Specification :

Antenna Trainer



Laboratory Instrument India