



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
FM Modulation and Demodulation Trainer

**Product Code :**  
LBNY-0002-610013



**Description :**

FM Modulation and Demodulation Trainer

**Technical Specification :**

With inbuilt Modulating Signal & Carrier signal.

Instrument comprises of fixed DC Regulated power supply of  $\pm 12$  Volts.

Built in modulating signal generator of 2KHz & 4KHz, Amplitude of 3V Peak to Peak.

Circuit Diagrams of modulator & demodulator printed on front panel.

2mm Sockets provided for connection on the front panel with standard accessories & instruction manual.

FM Modulation and Demodulation Trainer, FM Modulation and Demodulation Trainer Equipments, FM Modulation and Demodulation Trainer Tools, FM Modulation and Demodulation Trainer Tool Kits, FM Modulation and Demodulation Trainer Manufacturers, FM Modulation and Demodulation Trainer Suppliers from India, China, Kenya



**Laboratory Instrument India**

