



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Surveying Learning System

**Product Code :**  
LBNY-0005-1940004



**Description :**

The Surveying Learning System introduces the fundamentals of surveying, focusing on determining the terrestrial or three-dimensional position of points and the distances and angles between them; this signature approach to curriculum reinforces both theory and practice, which produces a well-rounded understanding of the topic.

These components are all industrial quality not only to ensure durability, but also to help learners become better prepared for what they will encounter on the job. All products are made from top-notch materials and carefully crafted (hand-welded, painted, and silk-screened) to create tough, attractive, well designed learning systems that facilitate learning and will serve teachers and students for years. The Surveying Learning System includes Level-Transit, optical type with 20x magnification, horizontal and vertical adjustments and locks, Level-Transit tripod with adjustable legs and mounting platform; Measuring Tape, fiberglass, 200-ft., (100) Surveying Flags, plastic film on wire stakes.

Because surveying is utilized in so many aspects of industry, the skills taught by the Surveying Learning System opens the doors to careers in various fields, such as manufacturing, transportation, and construction. After completing this training system, learners will be able to operate surveying and mapping tools, including global positioning systems (GPS) associated with surveying, but will also understand the relevant key concepts.

**Features:**

Introduces the fundamentals of surveying, focusing on determining the terrestrial or three-dimensional position of points and the distances and angles between them.

**Technical Specification :**

Surveying Learning System



**Laboratory Instrument India**