



Email : sales@laboratoryinstrumentindia.com

Product Name :
Stall Bar Learning System

Product Code :
LBNY-0005-1940008



Description :

The Stall Bar Learning System includes a pneumatic 50NM nutrunner with stall bar, 19MM socket, CD-ROM based student learning materials for both theory and lab, and teachers assessment guide. The Stall Bar Learning System teaches learners how to use a stall bar torque tool to assemble threaded fasteners. These are often used in assembly operations requiring very accurate torque. Stall Bar Learning System uses industrial quality components to help assure that students are better prepared for what they will encounter on the job and to withstand frequent student use. The CD-ROM based curriculum includes virtual simulation of key skills, allowing learners to practice skills in advance of actually working on the equipment. The CD-ROM based curriculum can be accessed locally or web delivered, allowing businesses to increase delivery of high quality assembly practices in far flung organizations.

Features:

- Industrial quality components
- Hands-on skills
- Teaches learners how to use a stall bar torque tool to assemble threaded fasteners.

Technical Specification :

Stall Bar Learning System

Laboratory Instrument India

