



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Oven 50 Liters Natural Convection

**Product Code :**  
LBNY-0005-1950003



#### **Description :**

The oven is equipped of a dual safety thermostat with higher thermic threshold to prevent accidental over-temperatures, and to ensure safe working conditions. Designed for drying, baking, conditioning and moisture determination. Sturdy manufacture, double walled with 60 mm thick glass fiber for thermal insulation. Natural convection, digital thermoregulator-indicator. Exterior front part is stainless steel made; while interior chamber, grid shelves and external walls are made from zinc coated steel. The oven is supplied complete with two grid shelves easily removable and that can be positioned at various heights, pilot light, and exhaust holes for fast cooling. Temperature from ambient to 200°C. is controlled by a digital thermoregulator-indicator.

#### **Technical Specification :**

Inside dimensions: 350x360x390 mm  
Outside dimensions: 590x460x620 mm  
Wattage: 750W  
Doors: n.1  
Power supply: 230V 50-60Hz 1ph  
Weight: 34 kg.



**Laboratory Instrument India**

