



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Oven 220 Liters Forced Ventilation Digital Thermostat

**Product Code :**  
LBNY-0005-1950004



#### **Description :**

##### **Main Features:**

Air outlet control by valves operated manually  
Regulation of fan speed (0-100%)  
High quality thermoinsulation material  
Forced ventilation airflow  
300°C maximum temperature  
Double over heat protection system  
Stainless steel chamber and trays with Silicone gasket  
Low power consumption  
PID electronic regulator, double digital display  
High temperature (uniformity  $\pm 2\%$  and precision  $\pm 0.3^\circ\text{C}$ ) fully compliant with Standards.

##### **Technical Specification :**

Doors: n.1  
Wattage: 4 kW  
Inside dimensions: 730x500x620 mm  
Outside dimensions: 930x880x915 mm  
Power supply: 230V 50-60Hz 1ph (3ph for yy and zz)  
Weight: 102 kg.



**Laboratory Instrument India**