



Email : sales@laboratoryinstrumentindia.com

Product Name :
Oven 220 Liters Forced Ventilation Digital Thermostat

Product Code :
LBNY-0005-1950004



Description :

Main Features:

Air outlet control by valves operated manually
Regulation of fan speed (0-100%)
High quality thermoinsulation material
Forced ventilation airflow
300°C maximum temperature
Double over heat protection system
Stainless steel chamber and trays with Silicone gasket
Low power consumption
PID electronic regulator, double digital display
High temperature (uniformity $\pm 2\%$ and precision $\pm 0.3^\circ\text{C}$) fully compliant with Standards.

Technical Specification :

Doors: n.1
Wattage: 4 kW
Inside dimensions: 730x500x620 mm
Outside dimensions: 930x880x915 mm
Power supply: 230V 50-60Hz 1ph (3ph for yy and zz)
Weight: 102 kg.



Laboratory Instrument India