



Email : sales@laboratoryinstrumentindia.com

Product Name :
Oven 420 Liters Forced Ventilation Digital Thermostat

Product Code :
LBNY-0005-1950005



Description :

Main Features:

Forced ventilation airflow
Regulation of fan speed (0-100%)
Stainless steel chamber and trays with Silicone gasket
PID electronic regulator, double digital display.
High quality thermo insulation material
Double over heat protection system
Low power consumption
300°C maximum temperature
High temperature (uniformity $\pm 2\%$ and precision $\pm 0.3^\circ\text{C}$) fully compliant with Standards
Air outlet control by valves operated manually

Technical Specification :

Doors: n.2
Wattage: 6.2 kW
Power supply: 230V 50-60Hz 1ph (3ph for yy and zz)
Inside dimensions: 1000x500x620 mm
Outside dimensions: 1200x930x1200 mm
Weight: 155 kg.



Laboratory Instrument India