




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Oven 100 Liters Natural Convection

**Product Code :**  
LBNY-0005-1950006



**Description :**

The oven is equipped of a dual safety thermostat with higher thermic threshold to prevent accidental over-temperatures, and to ensure safe working conditions. The oven is supplied complete with two grid shelves easily removable and that can be positioned at various heights, pilot light, and exhaust holes for fast cooling. Natural convection, digital thermoregulator-indicator. Designed for drying, baking, conditioning and moisture determination. Sturdy manufacture, double walled with 60 mm thick glass fiber for thermal insulation. Exterior front part is stainless steel made; while interior chamber, grid shelves and external walls are made from zinc coated steel. Temperature from ambient to 200°C. is controlled by a digital thermoregulator-indicator.

**Technical Specification :**

Doors: n.1  
Wattage: 1200W  
Power supply: 230V 50-60Hz 1ph  
Inside dimensions: 400x420x600 mm  
Outside dimensions: 640x515x805 mm  
Weight: 40 kg.



**Laboratory Instrument India**

