



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bench Top Anaerobic Digester

Product Code :
LBNY-0005-1960006



www.laboratoryinstrumentindia.com

Description :

These processes allow to reduce the sludge production and to produce fuel gas. The Bench Top Anaerobic Digester is designed as a comprehensive study facility of water treatment process involving bacteria in absence of air. The unit is supplied with manuals that describe each component of the system, installation and utilization procedures and give many exercises.

Technical Specification :

Required services:

Electrical power supply: 110/220V, single-phase, 50/60 Hz; 1kW.

Weight and dimensions:

Dimensions: 550x1000x800 h mm.

Weight (empty): 100kg.



Laboratory Instrument India