



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Still Plant

Product Code :
LBNY-0005-1960007



www.laboratoryinstrumentindia.com

Description :

The distillation tank is completely made of borosilicate glass with a low coefficient of thermal expansion. The electric heating requires a power of 4 kW for a production of approx 6 l/h with a conductivity of 1-2 $\mu\text{S/cm}$. Water Still Plant can be used both to produce distilled water for laboratory use and as a pilot plant for educational purposes. The plant is fitted with a safety alarm against malfunction or lack of water pressure. A training handbook completes the supply.

Experiments:

Study of the operation of a water condenser.
Distillation of a known quantity of water.
Effect of the water quantity variation.

Technical Specification :

Water Still Plant



Laboratory Instrument India

