



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Water Still Plant

**Product Code :**  
LBNY-0005-1960007



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

The distillation tank is completely made of borosilicate glass with a low coefficient of thermal expansion. The electric heating requires a power of 4 kW for a production of approx 6 l/h with a conductivity of 1-2  $\mu\text{S/cm}$ . Water Still Plant can be used both to produce distilled water for laboratory use and as a pilot plant for educational purposes. The plant is fitted with a safety alarm against malfunction or lack of water pressure. A training handbook completes the supply.

**Experiments:**

Study of the operation of a water condenser.  
Distillation of a known quantity of water.  
Effect of the water quantity variation.

**Technical Specification :**

Water Still Plant



**Laboratory Instrument India**

