



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bench Top Aerobic Digester

Product Code :
LBNY-0005-1960008



www.laboratoryinstrumentindia.com

Description :

A synthetically prepared waste water may be used to gain a working knowledge of the operational parameters and purification efficiencies. The Bench Top Aerobic Digester is designed as a comprehensive study facility of biological water treatment process by means of the continuous activated sludge. The unit is supplied with manuals that describe each component of the system, installation and utilization procedures and give many exercises.

Technical Specification :

Required services:

Electrical power supply: 110/220V, single-phase, 50/60 Hz; 1kW.

Weight and dimensions:

Dimensions: 550x1000x600 h mm

Weight (empty): 40kg.



Laboratory Instrument India

