




Email : sales@laboratoryinstrumentindia.com

Product Name :
Base Unit

Product Code :
LBNY-0005-19800012



Description :

The base unit can be set up horizontally, vertically and in inclined positions to suit experiment. The Base Unit is easy to set up with no assembly tools needed. Screen-printed design includes a measuring scale to ensure repeatable exercises.

Features:

- Can be assembled in vertical inclined & horizontal configurations
- Tool-less assembly
- Screen-printed design also includes a measuring scale to help identify correct length thumbscrews etc
- Sturdy common base for assembling all experiments to
- Improve the student's dexterity by self-assembly with graphical instructions provided
- Base Unit backboard has a screen-printed grid reference which all assembly instructions refer to for the location of parts
- Measuring scale included on backboard
- Flat-pack for easy storage
- Simple grid reference system to assist with assembly of experiments.

Technical Specification :

Base Unit

Laboratory Instrument India

