



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Aircraft Propeller Trainer

**Product Code :**  
LBNY-0005-1990004



**Description :**

Aircraft Propeller Trainer comprising a three blade propeller mounted on a free standing fixture. The range of propeller blade movement is demonstrated by the application of air pressure. The propellers are mounted on a moveable stand. The trainer equipped with a constant speed, feathering and reverse pitch propeller.

**Technical Specification :**

**Optional Item:**

The CBT may be used by students for self paced learning, and it may also be used by instructors to present high quality graphics and animations during the classroom lecture. The propeller trainer can be supplied with a high quality computer based training (CBT) course ware at an additional cost. The CBT covers the full scope of topics relevant to propellers, as taught in a typical aircraft maintenance training program, in detail with rich multimedia contents.



**Laboratory Instrument India**