



Email : sales@laboratoryinstrumentindia.com

Product Name :
Aircraft Turboprop Engine Trainer

Product Code :
LBNY-0005-1990008



Description :

The engine is operated from a semi enclosed test cell cabin that has all the necessary instrumentation and controls to run the engine. The engines are mounted with all necessary electrical wiring, engine accessories, controls, and instruments that are needed for running the engine. The trainer comes complete with an operational turboprop engine. The engine is connected to propeller assembly including spinner and governor. The engine trainer may be mounted on wheels for easy mobility. For fixed installations, the trainer is provided with heavy duty vibration and shear mounts to absorb vibrations. The students may also be trained to remove and replace parts for field level maintenance practice.

Technical Specification :

Exhaust:Exhaust Outlets

Propeller: Hartzell3 blade constant speed full feathering reversing propeller

Engine: Choices Pratt & Whitney Standard. Garrett or Walter optional

Electrical:Complete engine electrical system including Starter system, ignition leads,ignition unit, plugs, and battery.



Laboratory Instrument India

