



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Thrust Reverser System Trainer

**Product Code :**  
LBNY-0005-19900010



**Description :**

Thrust Reverser System Trainer is a functional system that represents the systems installed in turbine powered aircraft. The thrust reverser is pneumatically operated using the shop air above 50 psi. The operation is controlled by a throttle. It can be used by the instructors to demonstrate the actual operation of a thrust reverser, besides allowing hands-on operation like component replacement.

**Technical Specification :**

The trainer includes the following:  
Complete Thrust Reverser System  
System Arming Indicator  
MLG Air/Ground Sensor  
Control Throttle with reverse thrust gate  
System Control Panel  
Oil Pressure Sensor  
Fire Warning Switch  
Durable, sturdy steel frame with swivel caster wheels & brakes  
Safety Cage  
Power Supply (requires 110V AC or 220V AC input).

---

# Laboratory Instrument India

