



Email : sales@laboratoryinstrumentindia.com

Product Name :
Thrust Reverser System Trainer

Product Code :
LBNY-0005-19900010



Description :

Thrust Reverser System Trainer is a functional system that represents the systems installed in turbine powered aircraft. The thrust reverser is pneumatically operated using the shop air above 50 psi. The operation is controlled by a throttle. It can be used by the instructors to demonstrate the actual operation of a thrust reverser, besides allowing hands-on operation like component replacement.

Technical Specification :

The trainer includes the following:

- Complete Thrust Reverser System
- System Arming Indicator
- MLG Air/Ground Sensor
- Control Throttle with reverse thrust gate
- System Control Panel
- Oil Pressure Sensor
- Fire Warning Switch
- Durable, sturdy steel frame with swivel caster wheels & brakes
- Safety Cage
- Power Supply (requires 110V AC or 220V AC input).

Laboratory Instrument India

