



Email : sales@laboratoryinstrumentindia.com

Product Name :
Teardown Turbojet Engine

Product Code :
LBNY-0005-19900014



Description :

The tear down engine is used for hands-on training for the dis assembly and re-assembly of engine parts and sections. The engine stand allows easy positioning of the engine during work. The stand includes casters that allow convenience in moving the engine as desired. The students can also perform engine inspection. In conjunction with our engine cutaway, students can gain a thorough familiarization with engine construction and assembly.

Engine Features:

Length: 45 inches
Diameter: 17.7 inches
Type: Turbojet engine
Combustors: annular
Turbine: 2 stages
Dry weight: 410 Lbs. approx.
Compressor: 8 stages axial.

Technical Specification :

Teardown Turbojet Engine

Laboratory Instrument India

