



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutting Of Turbofan Engine Model

Product Code :
LBNY-0005-19900017



Description :

The Turbofan Engine Cutaway Model is powered by an AC motor to turn freely to demonstrate the phased operation of compressor and turbine.

Engine Features:

Fan/Compressor Stages: 1/8
Length (Inches): 75.5
Core / Power Turbine Stages: 2/1
Max Diameter (Inches): 33
Bypass Ratio: 2.0
Dry Weight (Lbs.): 725 - 737.

Cutaway Features:

Motorized cutaway on engine stand with casters
Cutaway provides a clear view of the:
Compressor Section
Combustion section
Turbine section.

Technical Specification :

Cutting Of Turbofan Engine Model



Laboratory Instrument India