



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Basic Electrical Trainer

**Product Code :**  
LBNY-0005-2000006



**Description :**

**Components:**

DC Power Switch with indicator light  
Master Power Switch with indicator light  
Solder-tail aircraft connector  
Audio alarm  
Circuit Breakers  
DC power buss  
Relays  
Switches  
Terminal strips  
Voltmeter  
Ammeter  
Indicator lamps  
Power Requirement: 110V AC or 220V AC  
Dimensions: 18 L x 18 W x 10 H inches.

**Technical Specification :**

Basic Electrical Trainer

---

# Laboratory Instrument India

