



Email : sales@laboratoryinstrumentindia.com

Product Name :
Basic Electrical Trainer

Product Code :
LBNY-0005-2000006



Description :

Components:

- DC Power Switch with indicator light
- Master Power Switch with indicator light
- Solder-tail aircraft connector
- Audio alarm
- Circuit Breakers
- DC power buss
- Relays
- Switches
- Terminal strips
- Voltmeter
- Ammeter
- Indicator lamps
- Power Requirement: 110V AC or 220V AC
- Dimensions: 18 L x 18 W x 10 H inches.

Technical Specification :

Basic Electrical Trainer

Laboratory Instrument India

