




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hydraulic Landing Gear Trainer

**Product Code :**  
LBNY-0005-2000008



#### **Description :**

The Landing Gear Mechanism Trainer includes all assemblies and components of the standard hydraulically operated landing gear system. Two types of trainers are available; comprises of three wheels comprising two main landing gears and one nose gear, and has one wheel. The Landing Gear Mechanism Trainer is very valuable equipment for the hands-on training for aircraft maintenance training schools. Trainer is a very effective hands on training system for landing gear system. Both models depict the landing gear mechanism as found in actual aircraft. The model is mounted on a move able stand to facilitate easy positioning of the trainer for a clear view from all directions. It is essential to ensure that trainees not only understand the theory but they are also well prepared to maintain the landing gear system in the best possible manner, and in compliance with the prevailing safety standards.

#### **Features:**

- 1-Tire Assembly
- Filter
- Fluid Reservoir
- Actuating Cylinder
- Wheel
- Tires
- 2-Hydraulic Brake Components
- Fluid Reservoir
- Brake Assembly
- Master Brake Cylinder
- 3- Hydraulically Operated Landing Gear
- Electrically Driven Hydraulic Pump

---

Relief Valves  
Emergency Hand Pump.

**Technical Specification :**

Hydraulic Landing Gear Trainer



**Laboratory Instrument India**