



Email : sales@laboratoryinstrumentindia.com

Product Name :
Antiskid Brake System Trainer

Product Code :
LBNY-0005-20000010



Description :

The Anti skid Brake System Trainer includes standard hydraulic brake system plus the anti skid assemblies and components. The system is designed to meet or exceed the FAR requirements for instruction on aircraft anti-skid braking system. The Anti skid braking system is an essential part of a modern jet aircraft. There is a definite need for training aircraft maintenance technicians in this area. Our Trainer is an effective hands-on training system for anti skid braking system. Hands-on training is the only effective method of training to ensure that trainees not only understand the theory but they are also well prepared to maintain the anti skid braking system in the best possible manner.

Technical Specification :

The anti skid brake system trainer consists of:

- Electronic Control Box
- Antiskid Control Valve
- Landing Gearwheel and Tire Assembly
- Skid Sensor
- Control Panel
- Motor-powered Hydraulic System.

Power Supply:

The trainer requires 110V AC 60 Hz.



Laboratory Instrument India