




Email : sales@laboratoryinstrumentindia.com

Product Name :
Aircraft Electrical Landing Gear Trainer

Product Code :
LBNY-0005-20000013



Description :

This display may be used for demonstration and for hands-on training assignments for the students. The nose gear is mounted on a sturdy steel frame, allowing full access to the landing gear.

Features:

1-Tire Assembly
Shimmy dampener
Nose steering mechanism
Shock absorber
Taxi Light
Oleo strut
Dimensions: 36 L x 32 D x 50 H inches

Technical Specification :

Manuals & Documentation:
Students manual
Technical reference manual
Operation/Maintenance manual
Instructors manual.

Laboratory Instrument India

