




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Aircraft Electrical Landing Gear Trainer

**Product Code :**  
LBNY-0005-20000013



**Description :**

This display may be used for demonstration and for hands-on training assignments for the students. The nose gear is mounted on a sturdy steel frame, allowing full access to the landing gear.

**Features:**

1-Tire Assembly  
Shimmy dampener  
Nose steering mechanism  
Shock absorber  
Taxi Light  
Oleo strut  
Dimensions: 36 L x 32 D x 50 H inches

**Technical Specification :**

Manuals & Documentation:  
Students manual  
Technical reference manual  
Operation/Maintenance manual  
Instructors manual.

---

# Laboratory Instrument India

