




Email : sales@laboratoryinstrumentindia.com

Product Name :
Cockpit Instrumentation Trainer

Product Code :
LBNY-0005-2010005



Description :

Three degrees of freedom instrument panel permits full demonstration of attitude and directional gyro functions. The Cockpit Instrumentation System with Electronic Flight Instruments and electronic engine display.

Standard Features:

- Functional engine monitoring system connected to engine sensors
- Provision of engine sensor simulation
- Functional digital primary flight display mounted on a panel that can simulate roll, pitch, and yaw movements controlled by a mechanism operated by a control yoke.
- Functional pitot-static system with provision for conducting pitot static system checks
- Standby Instruments installed.

Technical Specification :

The Model is fitted with the following sensors:

- EGT Thermo couples (Qty. 4)
- CHT Thermo couples (Qty. 4)
- Oil Temperature Sensor
- OAT Sensor
- Pitot-static Tube
- Fuel Flow Sensor
- Manifold Pressure Sensor
- Fuel Pressure Sensor
- Engine RPM Sensor

Fuel Level Sensor
Ammeter Shunt
Oil Pressure Sensor.



Laboratory Instrument India