



Email : sales@laboratoryinstrumentindia.com

Product Name :
Global Positioning Systems (GPS)

Product Code :
LBNY-0001-810005



Description :

Global Positioning Systems (GPS)

Technical Specification :

Global Positioning Systems (GPS)

The GPS device has a 2.2-inch or larger sunlight-readable color display with a resolution of 240 x 320 pixels. Preloaded with Active maps.

It will have support for GPS satellite systems and 8GB of internal memory for downloading maps, as well as a micro SD card slot.

Battery life up to 25 hours in GPS mode with standard size batteries.



Laboratory Instrument India