



Email : sales@laboratoryinstrumentindia.com

Product Name :
Aircraft DME System Trainer

Product Code :
LBNY-0005-20100012



Description :

Aircraft Distance Measuring Equipment (DME) System Trainer is designed to provide hands-on training on a real and operational equipment approved for aircraft installation. The trainer includes aircraft quality electrical shielded wiring, cables, connectors, and harness assemblies and DME antenna with coaxial cable connections. The equipment is mounted on a mobile workstation containing air craft DME unit with built-in indicator, and DME antenna, and a test set. The trainer includes a test panel console for bench testing of DME equipment.

Features:

Durable, steel frame, Vertical standing panel with sturdy casters
Built-in DCP over Supply
Simulation of Distance, and Speed
Avionics panel complete with Honeywell DME with built-in Indicator for distance, speed, and time
Accessories storage drawers & cabinet.

Technical Specification :

DME Trainer Specifications:

Deliver ables include the following items:
Power supply 28V, 10A Test Set for Transponder and DME
DME Bench Test panel with current meter, volt meter, and fused input power
CBT Admin Module with Database Connectivity.
Computer-based Training.
Test harness assembly
Test Set coax cables

Panel-mounted
Cables
DME antenna
Test Antenna
DME Fault Simulation Panel.



Laboratory Instrument India