



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Aircraft DME System Trainer

**Product Code :**  
LBNY-0005-20100012



#### **Description :**

Aircraft Distance Measuring Equipment (DME) System Trainer is designed to provide hands-on training on a real and operational equipment approved for aircraft installation. The trainer includes aircraft quality electrical shielded wiring, cables, connectors, and harness assemblies and DME antenna with coaxial cable connections. The equipment is mounted on a mobile workstation containing air craft DME unit with built-in indicator, and DME antenna, and a test set. The trainer includes a test panel console for bench testing of DME equipment.

#### **Features:**

Durable, steel frame, Vertical standing panel with sturdy casters  
Built-in DCP over Supply  
Simulation of Distance, and Speed  
Avionics panel complete with Honeywell DME with built-in Indicator for distance, speed, and time  
Accessories storage drawers & cabinet.

#### **Technical Specification :**

##### **DME Trainer Specifications:**

Deliver ables include the following items:  
Power supply 28V, 10A Test Set for Transponder and DME  
DME Bench Test panel with current meter, volt meter, and fused input power  
CBT Admin Module with Database Connectivity.  
Computer-based Training.  
Test harness assembly  
Test Set coax cables

---

Panel-mounted  
Cables  
DME antenna  
Test Antenna  
DME Fault Simulation Panel.



**Laboratory Instrument India**