




Email : sales@laboratoryinstrumentindia.com

Product Name :
Aircraft Autopilot System Trainer

Product Code :
LBNY-0005-20100016



Description :

It has a modular design. It may be installed and operated as a roll axis only autopilot, roll and pitch axis, or a roll, pitch, and yaw axis autopilot. Autopilot Trainer is based on a Century 2000 Autopilot. The advanced design provides very smooth and precise corrections in aircraft attitude. The flexible design provides the economy of expand ability. The complete Century 2000 is a true 3-axes, flight director/autopilot.

System Components:

Roll Axis Control

The complete roll axis autopilot features continuous roll/NAV/voltage monitors, heading select, heading hold, NAV, /LOC/BC anticipative intercept and track automatic soft track, about 45 degree intercept and all angle intercept with HSI.

System Components:

Directional Gyro/Indicator

Roll Servo

Autopilot Computer/Control Unit

Attitude Indicator / Flight Director

Moveable Control Surfaces

Turn & Slip Indicator

Altitude Pre-selector

Autopilot Annunciator

Rate of Climb Indicator

Altitude Sensor

Yaw Servo

Pitch Servo.

Technical Specification :**Power Requirements:**

Operates from either 110V / 60 Hz or 220V / 50 Hz mains power.



Laboratory Instrument India