




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Perforated Plate With Adjustable Stem (Swell Plate)

**Product Code :**  
LBNY-0005-2020004



**Description :**

Contact end of the stem is easily locked in place with a knurled nut. Perforated 5-7/8" (149mm) dia. base with adjustable stem.

**Features:**

Easy to insert and remove  
Adjustable stem for locking plate.

**Technical Specification :**

Perforated Plate With Adjustable Stem (Swell Plate)



**Laboratory Instrument India**