



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ductility Testing Apparatus

Product Code :
LBNY-0005-2020007



Description :

The machine consists of a carriage moving over a lead screw. An electric motor driven reduction gear unit ensures smooth constant speed and continuous operation. Designed to test three specimens simultaneously.

Technical Specification :

It is equipped with an electric pump circulator and heater. The temperature is controlled thermostatically. The entire assembly is mounted with a stainless steel lined water bath completely encased in metal bound hardwood. Two rates of travel i.e. 50 mm/min and 10 mm/min are provided. Suitable for operation on 230 V, 50 Hz, Single Phase, A.C. supply.

Complete With:

Thermometer IP 38C, Range: 23°C to 27°C
Ductility Mould, with Base Plate 3 Nos.



Laboratory Instrument India