




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Test For Determination Of Solubility

**Product Code :**  
LBNY-0005-2020009



**Description :**

Present set-up is self-contained water re-circulating unit, provided with a Rubber ring for Gooch crucible, Gooch crucible, Filter discs, fiber glass, 25 ml dia. Pack of 100. Including: filter flask 500 ml capacity, funnel for Gooch crucible.

**Purpose**  
Used for the determination of solubility.

**Technical Specification :**

Test For Determination Of Solubility



**Laboratory Instrument India**