



Email : sales@laboratoryinstrumentindia.com

Product Name :
Seepage Apparatus

Product Code :
LBNY-0005-20200011



Description :

It consists of a self-contained bench. Glass fronted tank, aluminium rear wall, 6 pressure tappings, adjustable overflow, low and high range flow-meters with digital readout, sump tank and pump, dye injector, model set and instruction manual. a set of models. 230V, 50Hz. This bench is used for experimental studies of seepage and permeability problems in porous media.

Technical Specification :

Seepage Apparatus



Laboratory Instrument India