



Email : sales@laboratoryinstrumentindia.com

Product Name :
Gas Welding Device Moveable

Product Code :
LBNY-0004-820001



Description :

Gas Welding Device Moveable, Gas Welding Device Moveable Equipments, Gas Welding Device Moveable Tools, Gas Welding Device Moveable Tool Kits, Gas Welding Device Moveable Manufacturers, Gas Welding Device Moveable Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

High-capacity Oxygen and Acetylene regulators, with a full 2 7/12" diam'
Single stage regulators with reverse flow check valves, brass tool handle and brass valve body
Precision torch handle with dual check valves, cutting attachment and cutting tip (size No. 0)
UL-listed cutting torch and rosebud heating tip (size No. 6)
3 welding nozzles (size No. 0, No. 2, No. 4), all with a swaged tip and individual spiral mixer
Flint striker tool and goggles with a No. 5 shade
15' welding hose, tip cleaning tools
Perfect for an Amateur Metalworker or a Seasoned Professional
All-brass torch handle and check valves
Brass and steel cutting torch attachment

Brass and copper welding/heating tip
Material: Brass Torch, Poly Carrier
Hose Length: 15-Feet
Cutting Capacity: 1/2"
Welding Capacity: 3/16"
Oxygen Regulator: CGA 540
Acetylene Regulator: CGA510

Accessory

Cutting attachment
Cutting Nozzle
Welding Pipe 3Ea
Welding pipe joint
Regulator
Twin-Welding Hoses (10m)
Goggles
Tip Cleaner
Spark lighter
Spanner
Carry box
Welding Mask



Laboratory Instrument India