



Email : sales@laboratoryinstrumentindia.com

Product Name :
CBR Test Apparatus Motorized

Product Code :
LBNY-0005-20200014



Description :

Upper beam can be adjusted in height. Foreseen of fast approach device of the base plate and electric end of stroke switches of the load plate to save the machine from wrong manipulations. CBR loading machine motorized 50 KN ASTM version 1,27 mm/min Load is applied through a screw jack driven by an electric motor at a constant penetration rate of 1,27 mm/min achieved by a built in gear box and assured also under load.

Technical Specification :

Load ring 50 KN capacity
Dial gauge holder
Dial gauge 10 x 0,01 mm
Power supply: 230 V 1ph 50 Hz 750 W
Dimensions: 430x380x1180 mm
Weight: 98 Kg.



Laboratory Instrument India