




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Apparatus For Distillation Of Cup Back Asphalt

**Product Code :**  
LBNY-0005-20200015



**Description :**

Complete with one distillation flask of 500 ml and one 100 ml crow receiver but without thermometer. Electric hot plate heated with energy regulator. The apparatus consists of a heavy gauge shield, lagged with asbestos cover in two parts and wire gauge. The shield is supported on a platform with chimney arranged for fitting to stand.

**Purpose:**

Used to measure the amount of the most volatile constituents in cut-back asphaltic products.

**Technical Specification :**

The apparatus consists of distillation flask, condenser tube, Shield and flask support, Bunsen burner, Cylinder receiver, adapter, corrosion resistant shield, thermometer -2 to + 400 0C subd. 1oC.



**Laboratory Instrument India**