



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
WS Workstation

**Product Code :**  
LBNY-0005-2030008



**Description :**

Castors are fitted to the Workstation which enable it to be easily moved between locations. The workstation can be used as a demonstration platform at the front of the classroom, but also doubles up as a storage unit for the experiment trays and Base Units.

**Features And Benefits:**

Ability to store both the experiment trays and base units  
The workstation can be configured to store:- Up to 12 base units and a maximum of 18 trays (2 base units occupy the space of 1 tray)  
18 trays if storing no base units  
12 base units and 12 trays  
Castors allow the workstation and its contents to be easily moved between locations  
Flat top can be used as a demonstration platform.

**Technical Specification :**

WS Workstation



**Laboratory Instrument India**

