



Email : sales@laboratoryinstrumentindia.com

Product Name :
WS Workstation

Product Code :
LBNY-0005-2030008



Description :

Castors are fitted to the Workstation which enable it to be easily moved between locations. The workstation can be used as a demonstration platform at the front of the classroom, but also doubles up as a storage unit for the experiment trays and Base Units.

Features And Benefits:

Ability to store both the experiment trays and base units
The workstation can be configured to store:- Up to 12 base units and a maximum of 18 trays (2 base units occupy the space of 1 tray)
18 trays if storing no base units
12 base units and 12 trays
Castors allow the workstation and its contents to be easily moved between locations
Flat top can be used as a demonstration platform.

Technical Specification :

WS Workstation



Laboratory Instrument India

