



Email : sales@laboratoryinstrumentindia.com

Product Name :
CO2 Welding Machine

Product Code :
LBNY-0004-820002



Description :

CO2 Welding Machine

Technical Specification :

CO2 welding machine
Specifications
CO2| Argon /Arc combined use
Rated voltage (eff): 1PW220Y/50-60H2
Rated input (KVA): 7,8-B
Utilization (7d): 60
Output voltage (A): TIG, MMA(10-120)/CU1(10-30)
No-load voltage (V): 52
Load voltage (V): 26.5
Welding ability (mm): 0.8-1.2
Accessory
Welding hammer, Brush
TIG Torch (5-6)M
Holder wire (5-6)M
Earth Sun (5-6)M
Pressure gauge
Plasma torch (5-6)M
Co2 Welding Wire: 0.9mm, 20kg 20Ea
Welding Mask.



Laboratory Instrument India