



Email : sales@laboratoryinstrumentindia.com

Product Name :
CO2 Welding Machine

Product Code :
LBNY-0004-820002



Description :

CO2 Welding Machine

Technical Specification :

CO2 welding machine
Specifications
CO2| Argon /Arc combined use
Ratedvoltageeff): 1PW220Y/50-60H2
Rated input(KVA): 7,8-B
Utilization(7d: 60
Output voltage(A): TIG, MMA(10-120)/CU1(10-30)
No-load voltage(V); 52
Load voltage(V): 26.5
Welding ability(mm): 0.8-1.2
Accessory
Welding hammer, Brush
TIGTorch(5-6)M
Holder wire (5-6)M
Earth Sun(5-6)M
Pressure gauge
Plasma torch (5-6)M
Co2 WeldingWile: 0.9mm,20kg20Ea
Welding Mask.



Laboratory Instrument India