



Email : sales@laboratoryinstrumentindia.com

Product Name :
Extension And Compression Of Springs Apparatus

Product Code :
LBNY-0005-17600029



Description :

Wall mounted apparatus (for loading spiral springs in either tension or compression with direct gravitational loads) with accurately divided scale and vernier.
Supplied complete with four compression and four extension springs each approximately 18 cm long (7") long.
Mass Hanger and Masses for (Hanger 1, Masses 13 x 1 kg)

Experimental Capability:

Modulus of rigidity for spring material.
Introduction to Spring Rate
Determines load and deflection relationships for close coiled springs.

Technical Specification :

Extension And Compression Of Springs Apparatus



Laboratory Instrument India