



Email : sales@laboratoryinstrumentindia.com

Product Name :
Rubber Block Apparatus

Product Code :
LBNY-0005-17600030



Description :

The comparator shows the degree of deformation of the block. Rubber Resistance Apparatus Consists of a panel on which a block of Rubber as well as a console equipped with a Dial are mounted. A total shear load of 12 kg is permitted.

Experimental Ability:

Determination of the stiffness modulus for the rubber
Observation of the driving force and the Shear deformation on a rubber block.
Determination of the variation in the deviation of the applied load
Demonstration of the relationship between the Shear and shear deformation.

Technical Specification :

Functionalities:

Wall unit
Allows the study of the modulus of stiffness of rubber.



Laboratory Instrument India

