



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Rubber Block Apparatus

**Product Code :**  
LBNY-0005-17600030



**Description :**

The comparator shows the degree of deformation of the block. Rubber Resistance Apparatus Consists of a panel on which a block of Rubber as well as a console equipped with ab Dial are mounted. A total shear load of 12 kg is permitted.

**Experimental Ability:**

Determination of the stiffness modulus for the rubber  
Observation of the driving force and the Shear deformation on a rubber block.  
Determination of the variation in the deviation of the applied load  
Demonstration of the relationship between the Shear and shear deformation.

**Technical Specification :**

**Functionalities:**

Wall unit  
Allows the study of the modulus of stiffness of rubber.



**Laboratory Instrument India**

