



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flood Lighting Set

Product Code :
LBNY-0005-1820003



www.laboratoryinstrumentindia.com

Description :

Flood Lighting Set is lighting equipment offered consists of four sealed beam flood lamps each housed in an aluminium body. Flood Lighting Set is essential to have the experimental area within the frame well lit. Flood Lighting Set lamp units slide in a 3.6 metre long track supported 1 metre from the front edge of the MAGNUS frame with cantilever brackets.

A mains connector and 6 metres of cable complete the assembly.

Specification:

A 3.6 m lighting track supplied at 1 m from top edge of MAGNUS frame

To be supplied with a mains connector and 6 m of cable

Four flood lamps in fittings that slide in the track.

Technical Specification :

Flood Lighting Set



Laboratory Instrument India

